

FGR Mouse mAb

Catalog No: #63638

Package Size: #63638-1 50ul #63638-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FGR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human FGR expressed in E. Coli.
Target Name	FGR
Other Names	SRC2; c-fgr; c-src2; p55c-fgr; p58c-fgr
Accession No.	P09769
Calculated MW	56kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using FGR mouse mAb against HL60
(1) Raw264.7 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human
Spleen tissues using FGR mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.