

FLT1 Mouse mAb

Catalog No: #63642

Package Size: #63642-1 50ul #63642-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FLT1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant extracellular fragment of human FLT1 fused with hlgGfc tag expressed in HEK293 cells.
Target Name	FLT1
Other Names	FLT; VEGFR1
Accession No.	P17948
Calculated MW	100kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ELISA:1/10000

Images



Western blot analysis using FLT1 mouse mAb against extracellular domain of human FLT1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.