

FLT3 Mouse mAb

Catalog No: #63643

Package Size: #63643-1 50ul #63643-2 100ul


Description

Product Name	FLT3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FLT3 (AA: 930-991) expressed in E. Coli.
Target Name	FLT3
Other Names	FLK2; STK1; CD135; FLK-2
Accession No.	P36888
Calculated MW	130kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using FLT3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using FLT3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.