

FOXD3 Mouse mAb

Catalog No: #63646

Package Size: #63646-1 50ul #63646-2 100ul


Description

Product Name	FOXD3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human FOXD3 expressed in E. Coli.
Target Name	FOXD3
Other Names	AIS1; HFH2; Genesis; FOXD3
Accession No.	Q9UJU5
Calculated MW	50kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using FOXD3 mouse mAb against NTERA-2 (1), HUVE-12 (2), HEK293 (3), Hela (4), Jurkat (5), K562 (6), RAW264.7 (7), NIH/3T3 (8), and COS7 (9) cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues (left) and ovarian cancer tissues (right) using FOXD3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.