

FOXP2 Mouse mAb

Catalog No: #63648

Package Size: #63648-1 50ul #63648-2 100ul


Description

Product Name	FOXP2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FOXP2 (AA: 641-740) expressed in E. Coli.
Target Name	FOXP2
Other Names	SPCH1; CAGH44; TNRC10
Accession No.	O15409
Calculated MW	80kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using FOXP2 mouse mAb against HepG2 (1) cell lysate.



Flow cytometric analysis of HeLa cells using FOXP2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.