

FRK Mouse mAb

Catalog No: #63650

Package Size: #63650-1 50ul #63650-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FRK Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of FRK (aa2-300) expressed in E. Coli.
Target Name	FRK
Other Names	GTK; RAK; PTK5
Accession No.	P42685
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human Breast tissues using FRK mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.