

FSHB Mouse mAb

Catalog No: #63652

Package Size: #63652-1 50ul #63652-2 100ul


Description

Product Name	FSHB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FSHB (AA: 19-129) expressed in E. Coli.
Target Name	FSHB
Accession No.	P01225
Calculated MW	14.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using FSHB mouse mAb against Cos7 (1), HepG2 (2), HEK293 (3) cell lysate.



Flow cytometric analysis of HepG2 cells using FSHB mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.