

FSHB Mouse mAb

Catalog No: #63653

Package Size: #63653-1 50ul #63653-2 100ul


Description

Product Name	FSHB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FSHB (AA: 2x 19-129) expressed in E. Coli.
Target Name	FSHB
Other Names	HH24
Accession No.	P01225
Calculated MW	14.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using FSHB mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human hypophysis tissues using FSHB mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.