

FYN Mouse mAb

Catalog No: #63656

Package Size: #63656-1 50ul #63656-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FYN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Purified recombinant fragment of human FYN expressed in E. Coli.
Target Name	FYN
Other Names	SLK; SYN; MGC45350
Accession No.	P06241
Calculated MW	61kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using FYN mouse mAb against NIH/3T3 (1) and Hela (2) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.