

FZD5 Mouse mAb

Catalog No: #63657

Package Size: #63657-1 50ul #63657-2 100ul


Description

Product Name	FZD5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human FZD5 (AA:151-217) expressed in E. Coli.
Target Name	FZD5
Other Names	HFZ5; C2orf31; MGC129692; DKFZp434E2135
Accession No.	Q13467
Calculated MW	64.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using FZD5 mouse mAb against A549 (1), and PC-3 (2) cell lysate.



Flow cytometric analysis of HepG2 cells using FZD5 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.