

G6D Mouse mAb

Catalog No: #63658

Package Size: #63658-1 50ul #63658-2 100ul


Description

Product Name	G6D Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human G6D (AA: 20-104) expressed in E. Coli.
Target Name	G6D
Other Names	LY6G6D; NG25; LY6-D; MEGT1; C6orf23
Accession No.	O95868
Calculated MW	13.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using G6D mouse mAb against Jurkat (1), MCF-7 (2), A549 (3),K562 (4),Hela (5),C6 (6),COS-7 (7),and NIH/3T3 (8) cell lysate.



Immunofluorescence analysis of HeLa cells using G6D mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Jurkat cells using G6D mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.