

GAB1 Mouse mAb

Catalog No: #63660

Package Size: #63660-1 50ul #63660-2 100ul


Description

Product Name	GAB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GAB1 expressed in E. Coli.
Target Name	GAB1
Accession No.	Q13480
Calculated MW	76.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);


Western blot analysis using GAB1 mouse mAb against HEK293 (1) cell lysate.



Flow cytometric analysis of Jurkat cells using GAB1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using GAB1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human kidney tissues using GAB1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.