

GABBR2 Mouse mAb

Catalog No: #63662

Package Size: #63662-1 50ul #63662-2 100ul


Description

Product Name	GABBR2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GABBR2 (AA: 319-483) expressed in E. Coli.
Target Name	GABBR2
Other Names	HG20; GPR51; GPRC3B; GABABR2; HRIHFB2099
Accession No.	O75899
Calculated MW	105.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using GABBR2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Secondary antibody from Fisher (Cat#: 35503)



Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using GABBR2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.