

GABPA Mouse mAb

Catalog No: #63663

Package Size: #63663-1 50ul #63663-2 100ul


Description

Product Name	GABPA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human GABPA (aa120-190) expressed in E. Coli.
Target Name	GABPA
Other Names	NFT2; NRF2
Accession No.	Q06546
Calculated MW	51kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using GABPA mouse mAb against Hela (1), A549 (2), MCF-7 (3), NIH/3T3 (4) and SMMC-7721 (5) cell lysate.



Confocal Immunofluorescence analysis of Hela cells using GABPA mouse mAb (green). Red: Actin filaments have been labeled using DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.