

GATA1 Mouse mAb

Catalog No: #63665

Package Size: #63665-1 50ul #63665-2 100ul


Description

Product Name	GATA1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GATA1 expressed in E. Coli.
Target Name	GATA1
Other Names	GF1; GF-1; NFE1; XLTT; ERYF1; GATA-1; GATA1
Accession No.	P15976
Calculated MW	42.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using GATA1 mouse mAb against K562 (1) cell lysate.



Immunofluorescence analysis of K562(left) cells using GATA1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human pancreatic cancer, using GATA1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.