

GATA3 Mouse mAb

Catalog No: #63666

Package Size: #63666-1 50ul #63666-2 100ul


Description

Product Name	GATA3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GATA3 (AA: (294-443)) expressed in E. Coli.
Target Name	GATA3
Other Names	HDR; HDRS
Accession No.	P23771
Calculated MW	48kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using GATA3 mouse mAb against T47D (1) cell lysate.



Flow cytometric analysis of Hela cells using GATA3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using GATA3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.