

GATA4 Mouse mAb

Catalog No: #63667

Package Size: #63667-1 50ul #63667-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GATA4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human GATA4 expressed in E. Coli.
Target Name	GATA4
Other Names	MGC126629
Accession No.	P43694
Calculated MW	46kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using GATA4 mouse mAb against rat fetal heart (1) and adult heart (2) tissues lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.