

Glucose-6-phosphate isomerase Mouse mAb

Catalog No: #63673

Package Size: #63673-1 50ul #63673-2 100ul


Description

Product Name	Glucose-6-phosphate isomerase Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GPI expressed in E. Coli.
Target Name	Glucose-6-phosphate isomerase
Other Names	AMF; NLK; PGI; PHI; GNPI; SA-36; GPI
Accession No.	P06744
Calculated MW	63kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 ELISA:1/10000


Images



Western blot analysis using GPI mouse mAb against HepG2 (1), SMMC-7721 (2) cell lysate and rat liver tissues lysate (3).



Confocal Immunofluorescence analysis of L-02 cells using GPI mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human Kidney tissues using GPI mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.