

GLUL Mouse mAb

Catalog No: #63675

Package Size: #63675-1 50ul #63675-2 100ul


Description

Product Name	GLUL Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human GLUL (AA: 2-121) expressed in E. Coli.
Target Name	GLUL
Other Names	GS; GLNS; PIG43; PIG59
Accession No.	P15104
Calculated MW	42kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using GLUL mouse mAb against Jurkat (1), and mouse liver (2) cell lysate.



Flow cytometric analysis of Hela cells using GLUL mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebellar tissues using GLUL mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GLUL mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.