

GMNN Mouse mAb

Catalog No: #63677

Package Size: #63677-1 50ul #63677-2 100ul


Description

Product Name	GMNN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human GMNN (AA: FULL 1-209) expressed in E. Coli.
Target Name	GMNN
Other Names	Gem; MGORS6
Accession No.	O75496
Calculated MW	23.6KDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of C6 cells using GMNN mouse mAb (green) and negative control (red).



Flow cytometric analysis of Hela cells using GMNN mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using GMNN mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using GMNN mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.