

GNAS Mouse mAb

Catalog No: #63679

Package Size: #63679-1 50ul #63679-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GNAS Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GNAS (AA: 42-188) expressed in E. Coli.
Target Name	GNAS
Other Names	AHO; GSA; GSP; POH; GPSA; NESP; GNAS1; PHP1A; PHP1B; PHP1C; C20orf45
Accession No.	P63092
Calculated MW	45.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A549 cells using GNAS mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using GNAS mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissues using GNAS mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.