

GOLGA2 Mouse mAb

Catalog No: #63681

Package Size: #63681-1 50ul #63681-2 100ul


Description

Product Name	GOLGA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GOLGA2 (AA: 1-205) expressed in HEK293-6e cells supernatant.
Target Name	GOLGA2
Other Names	GM130
Accession No.	Q08379
Calculated MW	113kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/50 - 1/200FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using GOLGA2 mouse mAb against HepG2 (1), Hela (2), K562 (3), and HEK293 (4) cell lysate.




Immunofluorescence analysis of HeLa cells using GOLGA2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using GOLGA2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using GOLGA2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using GOLGA2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.