

GPR19 Mouse mAb

Catalog No: #63685

Package Size: #63685-1 50ul #63685-2 100ul


Description

Product Name	GPR19 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GPR19 expressed in E. Coli.
Target Name	GPR19
Accession No.	Q15760
Calculated MW	47.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of C6 cells using GPR19 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.