

GRIA2 Mouse mAb

Catalog No: #63688

Package Size: #63688-1 50ul #63688-2 100ul


Description

Product Name	GRIA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GRIA2 (AA: 652-807) expressed in E. Coli.
Target Name	GRIA2
Other Names	GluA2; GluR2; gluR-B; GluR-K2
Accession No.	P42262
Calculated MW	98.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using GRIA2 mouse mAb against HeLa (1) cell lysate.



Flow cytometric analysis of SK-N-SH cells using GRIA2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.