

GRIK3 Mouse mAb

Catalog No: #63691

Package Size: #63691-1 50ul #63691-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GRIK3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GRIK3 (AA: extra 32-173) expressed in E. Coli.
Target Name	GRIK3
Other Names	EAA5; GLR7; GLUR7; GluK3; GluR7a
Accession No.	Q13003
Calculated MW	104kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using GRIK3 mouse mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunofluorescence analysis of HeLa cells using GRIK3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.