

GSK3 alpha Mouse mAb

Catalog No: #63702

Package Size: #63702-1 50ul #63702-2 100ul


Description

Product Name	GSK3 alpha Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of GSK3 alpha expressed in E. Coli.
Target Name	GSK3 alpha
Other Names	DKFZp686D0638; GSK3A
Accession No.	P49840
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung carcinoma (left) and skin carcinoma (right), showing cytoplasmic localization using GSK3 alpha mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.