

GSK3 alpha Mouse mAb

Catalog No: #63703

Package Size: #63703-1 50ul #63703-2 100ul


Description

Product Name	GSK3 alpha Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of GSK3 alpha expressed in E. Coli.
Target Name	GSK3 alpha
Other Names	DKFZp686D0638; GSK3A
Accession No.	P49840
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using GSK3 alpha mouse mAb against HeLa (1) and PC-3 cell lysate.



Immunofluorescence analysis of HeLa cells using GSK3 alpha mouse mAb showing cytoplasmic localization.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.