

GSN Mouse mAb

Catalog No: #63704

Package Size: #63704-1 50ul #63704-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GSN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GSN (AA: (673-783)) expressed in E. Coli.
Target Name	GSN
Other Names	ADF; AGEL
Accession No.	P06396
Calculated MW	85.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of Jurkat cells using GSN mouse mAb (green) and negative control (red).



Flow cytometric analysis of THP-1 cells using GSN mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using GSN mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectal cancer tissues using GSN mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.