

GST Tag Mouse mAb

Catalog No: #63705

Package Size: #63705-1 50ul #63705-2 100ul


Description

Product Name	GST Tag Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IP;IF
Species Reactivity	Tag
Immunogen Description	Purified fusion protein with GST protein tag expressed in E. Coli.
Target Name	GST Tag
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000IP:45677IHC:1/100 - 1/1000ICC:1/100 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images



Western blot analysis using GST mouse mAb against various fusion protein with GST tag.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.