

GYS1 Mouse mAb

Catalog No: #63710

Package Size: #63710-1 50ul #63710-2 100ul


Description

Product Name	GYS1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human GYS1 expressed in E. Coli.
Target Name	GYS1
Other Names	GSY; GYS; GYS1
Accession No.	P13807
Calculated MW	85kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using GYS1 mouse mAb against HeLa (1) and HEK293 (2) cell lysate.



Flow cytometric analysis of K562 cells using GYS1 mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.