

HAUSP Mouse mAb

Catalog No: #63716

Package Size: #63716-1 50ul #63716-2 100ul


Description

Product Name	HAUSP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HAUSP expressed in E. Coli.
Target Name	HAUSP
Other Names	TEF1; HAUSP; USP7
Accession No.	Q93009
Calculated MW	128kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using HAUSP mouse mAb against MCF-7 (1) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.