

HDAC10 Mouse mAb

Catalog No: #63718

Package Size: #63718-1 50ul #63718-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HDAC10 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HDAC10 (AA: 18-219) expressed in E. Coli.
Target Name	HDAC10
Other Names	HD10
Accession No.	Q969S8
Calculated MW	71.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using HDAC10 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using HDAC10 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.