

HDAC6 Mouse mAb

Catalog No: #63720

Package Size: #63720-1 50ul #63720-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HDAC6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HDAC6 (AA: 482-800) expressed in E. Coli.
Target Name	HDAC6
Other Names	HD6; JM21; CPBHM; PPP1R90
Accession No.	Q9UBN7
Calculated MW	131.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using HDAC6 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.