

HEXA Mouse mAb

Catalog No: #63722

Package Size: #63722-1 50ul #63722-2 100ul


Description

Product Name	HEXA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HEXA expressed in E. Coli.
Target Name	HEXA
Other Names	TSD
Accession No.	P06865
Calculated MW	60.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using HEXA mouse mAb against L1210 (1), and HL7702 (2) cell lysate.



Flow cytometric analysis of HeLa cells using HEXA mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.