

HFE Mouse mAb

Catalog No: #63723

Package Size: #63723-1 50ul #63723-2 100ul


Description

Product Name	HFE Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HFE expressed in E. Coli.
Target Name	HFE
Other Names	HH; HFE1; HLA-H; MVCD7; TFQTL2; MGC103790; dJ221C16.10.1
Accession No.	Q30201
Calculated MW	40kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunofluorescence analysis of HepG2 cells using HFE mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.