

HH3 Mouse mAb

Catalog No: #63724

Package Size: #63724-1 50ul #63724-2 100ul


Description

Product Name	HH3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Synthesized peptide fragment of human HH3 (AA: 121-136) expressed in E. Coli.
Target Name	HH3
Other Names	HIST2H3A;H3/n; H3/o
Accession No.	Q71DI3
Calculated MW	15.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using HH3 mouse mAb against K562 (1), C6(2),HEK293(3),PC-12(4) and NIH/3T3(5) cell lysate.



Flow cytometric analysis of NIH/3T3 cells using HH3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.