

His Tag Mouse mAb

Catalog No: #63727

Package Size: #63727-1 50ul #63727-2 100ul


Description

Product Name	His Tag Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IP;IHC;IF;FC
Species Reactivity	Tag
Immunogen Description	Synthesized peptide of His Tag (AA: HHHHHHHH-KLH).
Target Name	His Tag
Other Names	6 His tag; HHHHHH tag
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500-1/1000IP:1/20-1/200IHC:1/100 - 1/1000ICC:1/100 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images



Red: Control Antigen (100ng) Purple: Antigen (10ng), Green: Antigen (50ng) Blue: Antigen (100ng),



Western blot analysis using His tag mAb against recombinant protein. lane 1 :(100 ng); lane 2 :(50 ng); lane 3 :(25 ng); lane 4 :(10 ng); lane 5 :(2.5 ng); (Expected MW is 50 kDa).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.