

HIST2H3C(27Ac) Mouse mAb

Catalog No: #63728

Package Size: #63728-1 50ul #63728-2 100ul


Description

Product Name	HIST2H3C(27Ac) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Synthesized peptide of human HIST2H3C (AA: ATKAARK(Ac)SAPATGGV).
Target Name	HIST2H3C(27Ac)
Other Names	H3; H3.2; H3/M; H3F2; H3FM; H3FN
Accession No.	Q71DI3
Calculated MW	15.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using HIST2H3C(27Ac) mouse mAb against Hela (1), Lncap (2), Jurkat (3), and Jurkat (4) cell lysate.



Flow cytometric analysis of Raji cells using HIST2H3C(27Ac) mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using HIST2H3C(27Ac) mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using HIST2H3C(27Ac) mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.