

HLA-DRB1 Mouse mAb

Catalog No: #63738

Package Size: #63738-1 50ul #63738-2 100ul


Description

Product Name	HLA-DRB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HLA-DRB1 (AA: 30-266) expressed in E. Coli.
Target Name	HLA-DRB1
Other Names	SS1; DRB1; HLA-DRB; HLA-DR1B
Accession No.	P01911
Calculated MW	29.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using HLA-DRB1 mouse mAb against Raji (1) cell lysate.



Flow cytometric analysis of Raji cells using HLA-DRB1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human tonsil tissues using HLA-DRB1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.