

HLA-F Mouse mAb

Catalog No: #63740

Package Size: #63740-1 50ul #63740-2 100ul


Description

Product Name	HLA-F Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human HLA-F (AA: 22-305) expressed in E. Coli.
Target Name	HLA-F
Other Names	HLAF; CDA12; HLA-5.4; HLA-CDA12
Accession No.	P30511
Calculated MW	39kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using HLA-F mouse mAb against Mouse Liver (1) and Rat Liver (2) tissue lysate.



Flow cytometric analysis of THP-1 cells using HLA-F mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.