

HMOX1 Mouse mAb

Catalog No: #63741

Package Size: #63741-1 50ul #63741-2 100ul


Description

Product Name	HMOX1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HMOX1 (AA: 1-110) expressed in E. Coli.
Target Name	HMOX1
Other Names	HO-1; HSP32; HMOX1D; bK286B10
Accession No.	P09601
Calculated MW	32.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using HMOX1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human Liver tissues using HMOX1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using HMOX1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.