

HNF4A Mouse mAb

Catalog No: #63743

Package Size: #63743-1 50ul #63743-2 100ul


Description

Product Name	HNF4A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4552
Isotype	Mouse IgG2b
Applications	WB,IHC,IF,FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HNF4A (AA: 1-150) expressed in HEK293-6e cells supernatant.
Target Name	HNF4A
Other Names	TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7; HNF4a8; HNF4a9; NR2A21; TCF-14; HNF4alpha
Accession No.	P41235
Calculated MW	53 kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HepG2 cells using HNF4A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using HNF4A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using HNF4A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.