

HPRT1 Mouse mAb

Catalog No: #63746

Package Size: #63746-1 50ul #63746-2 100ul


Description

Product Name	HPRT1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HPRT1 (AA: FULL(1-218)) expressed in E. Coli.
Target Name	HPRT1
Other Names	HPRT; HGPRT
Accession No.	P00492
Calculated MW	24.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using HPRT1 mouse mAb against HeLa (1), A431 (2), A549 (3) cell lysate.



Flow cytometric analysis of Hela cells using HPRT1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human kidney tissues using HPRT1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.