

HTR3B Mouse mAb

Catalog No: #63756

Package Size: #63756-1 50ul #63756-2 100ul


Description

Product Name	HTR3B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human HTR3B (AA: extra 23-150) expressed in E. Coli.
Target Name	HTR3B
Other Names	5-HT3B
Accession No.	O95264
Calculated MW	50.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/100 - 1/500FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using HTR3B mouse mAb against SK-N-SH (1) and Hela (2) cell lysate.



Immunofluorescence analysis of HeLa cells using HTR3B mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using HTR3B mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.