

human IgG Mouse mAb

Catalog No: #63758

Package Size: #63758-1 50ul #63758-2 100ul


Description

Product Name	human IgG Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Human IgG was isolated from human sera and purified by chromatography.
Target Name	human IgG
Calculated MW	50kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using IgG mouse mAb (lane 1, 2) and Albumin mouse mAb (lane 3, 4) against human serum (lane 1, 3) and plasma (lane 2, 4).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.