

human IgG (Fc Specific) Mouse mAb

Catalog No: #63759

Package Size: #63759-1 50ul #63759-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	human IgG (Fc Specific) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified fusion protein with human IgG (Fc Specific) tag.
Target Name	human IgG (Fc Specific)
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using human IgG (Fc specific) mouse mAb against human serum (1).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.