

human IgG(Fc Specific) Mouse mAb

Catalog No: #63760

Package Size: #63760-1 50ul #63760-2 100ul


Description

Product Name	human IgG(Fc Specific) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Purified fusion protein with human IgG(Fc Specific) tag.
Target Name	human IgG(Fc Specific)
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using human IgG (Fc specific) mouse mAb against different fusion proteins with human IgG(Fc specific) tag.



Confocal Immunofluorescence analysis of HEK293 cells transfected with recombinant plasmid with human IgG Fc tag using human IgGFc mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.