

Human P16 Mouse mAb

Catalog No: #63761

Package Size: #63761-1 50ul #63761-2 100ul


Description

Product Name	Human P16 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of P16 expressed in E. Coli.
Target Name	Human P16
Other Names	Human P16
Accession No.	P42771
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung adenocarcinoma (A), esophageal squamous cell carcinoma (B), hepatic cell carcinoma (C), thyroid tumor (D), breast adenofibroma (E), breast infiltrating ductal carcinoma (F), normal cerebrum tissue (G), normal colon tissue (H), normal esophageal tissue (I), showing nuclear localization using P16 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human spleen tissues using P16 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.