

human Splunc2 Mouse mAb

Catalog No: #63762

Package Size: #63762-1 50ul #63762-2 100ul


Description

Product Name	human Splunc2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human Splunc2 expressed in E. Coli.
Target Name	human Splunc2
Other Names	BPIFA2; PSP; SPLUNC2; C20orf70; bA49G10.1
Accession No.	Q96DR5
Calculated MW	27kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of MCF-7 cells using human Splunc2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.