

ID2 Mouse mAb

Catalog No: #63763

Package Size: #63763-1 50ul #63763-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ID2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ID2 (AA: 1-134) expressed in E. Coli.
Target Name	ID2
Other Names	GIG8; ID2A; ID2H; bHLHb26
Accession No.	Q02363
Calculated MW	15kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of SK-N-SH cells using ID2 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human breast cancer tissues using ID2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ID2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.